

## Real-time PCR Diagnostic Kit for Rapid Identification of *Staphylococcus aureus*

Product No. MX-1601

### Intended Use

This product is the Real-time PCR Diagnostic Kit for Rapid Identification of *Staphylococcus aureus*, which can specifically amplify the specific DNA nucleic acid fragments of *Staphylococcus aureus* in the sample, and determine whether the sample contains *Staphylococcus aureus* by Ct value.

Specification 24 Reactions

### Principle and Interpretation

This kit uses real-time fluorescent PCR technology to perform specific amplification of *Staphylococcus aureus* by specific primers of *Staphylococcus aureus*. The fluorescence-labeled probes were used to label and track the amplification products, the reaction process was monitored online in real time, and whether the products were amplified or not could be analyzed in combination with corresponding software.

### Kit Contents

No.	Components	Amount
1	PCR Master mix (Powder)	24 test /Tube
2	RNase-Free ddH <sub>2</sub> O	1 mL/Tube
3	Lysis Buffer	5 mL/ Bottle
4	0.2 mL 8-Tube PCR Strips	3

### Self-provided materials

Sterilized 1.5mL centrifuge tube; sterilized 0.2mL PCR tube and tip; crushed ice or ice box; micropipette; centrifuge; vortex oscillator; metal bath; sterilized ultrapure water.

### Directions

#### 1. Sample Preparation:

Method: Reference to this manual:

a) Samples from culture: Take 1mL of the bacteria culture into a 1.5mL tube, centrifuge at 3000 ×g for 10min or 10000 ×g for 2min, discard the supernatant, add 200μL lysis buffer, mix well by vortex, heat at 99°C in metal bath or boil in a water bath for 10 minutes, cooling down in ice bath, centrifuge at 12000 ×g for 2min, and take supernatant to a new clean tube as DNA

template.

b) Samples from colony: Take one colony using the inoculating loop, suspend in 200 $\mu$ L lysis solution, mixed by vortex, heat at 99 $^{\circ}$ C in metal bath or boil in a water bath for 10 minutes, cooling down in ice bath, centrifuge at 12000 $\times$ g for 2 min, and take supernatant to a new clean tube as DNA template.

## 2. Prepare the PCR Reaction Mix

① Take out the PCR Master mix (Powder), add 575 $\mu$ L RNase-Free ddH<sub>2</sub>O to the tube, mix well by pipette or vortex.

Note: If more PCR mix remain, recommended aliquot to the 0.2mL PCR tube for storage for long time use.

② Take 23 $\mu$ L PCR mix to 0.2mL PCR tube.

③ Add 2 $\mu$ L DNA template from Step 1(Sample Preparation) into PCR tubes, with a total reaction volume of 25 $\mu$ L.

④ The corresponding positive control, negative control and sterile deionized water were added to the positive control, negative control and blank control systems respectively.

## 3. Reaction Condition

Stage	Cycles	Temperature	Time	Step	Fluorescent Signal# Acquisition
Pre-mutability	1	95 $^{\circ}$ C	5min	pre-mutability	No
Real-time PCR	40	95 $^{\circ}$ C	5sec	mutability	No
		60 $^{\circ}$ C	30-60 sec*	metallurgy	Yes

\*: When using different types of instruments for time setting, please follow the requirements of the instrument instruction manual for experimental operation, generally set at 30 sec.

#: The *Staphylococcus aureus*-specific gene has a fluorescent motif of FAM and a quenching motif of TAMRA.

## Result Analysis

Blank control: no FAM fluorescence signal was detected and no typical amplification curve appeared.

Negative control: no FAM fluorescence signal was detected and no typical amplification curve was observed.

Positive control (with causative gene): a FAM fluorescent signal was detected, a typical amplification curve was present, and the Ct value was <30.0.

All of the above need to be satisfied at the same time, otherwise this experiment is invalid.

Sample test results:

Ct values  $\geq 35$  were determined to be negative;

Samples were judged to be positive with a Ct value <30 and a typical amplification curve;

If the Ct value is  $30 \leq Ct < 35$  and a typical amplification curve is present, the DNA is reextracted for testing. If the result is still  $30 \leq Ct \text{ value} < 35$ , the sample is judged as positive; Ct value  $\geq$

35, the result is judged as negative.

### **Limitations of Method**

The target sequences detected by this kit are the conserved regions of specific genes in *Staphylococcus aureus*, which are highly conserved and stable. If the bacteria have genetic mutations at the target sequence, false negative results may occur. In addition, the procedure of sample collection, preparation, transportation and storage will also affect the test results.

### **Storage Conditions**

Storage at 2-8 °C, or store at -20 °C when not used for a long time.

### **Shelf Life**

24 month.

### **Announcements**

1. Please read the instruction of this kit carefully before the experiment and strictly follow the operation steps.
2. All components in this kit should be fully melted and mixed before use, and the high speed centrifugation for a short time is necessary.
3. This kit must be stored away from light, and the centrifuge tube and Tips with the DNase and RNase free should be autoclave before used. The whole operation process and the PCR laboratory should comply with the requirements of regulations such as “Administrative Measures for Clinical Gene Amplification Testing Laboratories in Medical Institutions” and “Working Guidelines for Clinical Gene Amplification Testing Laboratories in Medical Institutions” issued by the NHFPC. The waste and amplification products produced during the test should be properly treated to prevent cross-contamination.

*For laboratory use in industry or R&D purpose. Not for drug, household or other uses.*